

JPET #224246

**Supplemental Data:**

**Title:** Dual Inhibition of IL-23 and IL-17 Offers Superior Efficacy in Mouse Models of Autoimmunity

**Authors:** Paul R. Mangan, Linhui Julie Su, Victoria Jenny, Andrea L. Tatum, Caryn Picarillo, Stacey Skala, Noah Ditto, Zheng Lin, XiaoXia Yang, Pete Z. Cotter, David J. Shuster, Yunling Song, Virna Borowski, Rochelle L. Thomas, Elizabeth M. Heimrich, Brigitte Devaux, Ruchira Das Gupta, Irvith Carvajal, Kim W. McIntyre, Jenny Xie, Qihong Zhao, Mary Struthers, Luisa M. Salter-Cid

**Journal Name:** Journal of Pharmacology and Experimental Therapeutics

## **Supplemental Figure and Table Legends:**

### **Supplemental Table 1: Biochemical analysis of ATI-1221 and ATI-1249**

Kinetics for ATI-1221 and ATI-1249 (C-terminally PEGylated version of ATI-1221) were assessed against directly immobilized murine IL-23 using a Proteon instrument. Murine IL-23 at 2.5 to 5  $\mu\text{g/mL}$  in acetate pH 5.5 was immobilized on a GLC chip to a surface density of  $\sim 1000\text{-}2000$  RU.

Concentration series of adnectins in PBS-T were injected and dissociation monitored for up to 3 hours at 37C. Kinetic parameters were evaluated using a 1:1 fitting model with local Rmax. Data are the averaged kinetic parameters from the number of experiments indicated in the table.

### **Supplemental Figure 1: ATI-1221 and ATI-1249 Do Not Inhibit IL-12-induced Th1**

#### **Differentiation and Secretion of IFN- $\gamma$**

To differentiate Th1 cells, naive CD4<sup>+</sup>CD62L<sup>+</sup> T cells were enriched from the pooled spleen and lymph nodes of C57/B6 mice (Charles River Laboratories) by magnetic bead selection (Miltenyi Biotec) and seeded on 96 well plates coated with anti-CD3 (145-2C11, generated at BMS) in presence of 10 ng/ml recombinant mouse IL-12 (R&D Systems), 1  $\mu\text{g/ml}$  anti-CD28 antibody (3N7, generated at BMS) and 2  $\mu\text{g/ml}$  of anti-IFN- $\gamma$  (XMG1.2, eBioscience) neutralizing antibody. To test the activity of ATI-1221 and ATI-1249 for inhibiting IL-12-dependent responses, a range of each inhibitor as well as mouse anti-IL-12p40 antibody (R&D Systems) was added to the cultures during Th1 differentiation. After 4 days, the conditioned media from the cultures was harvested and assayed for IFN- $\gamma$  concentrations by ELISA (R&D Systems). The IL-12-dependent response was measured by calculating the difference between the level of cytokine induced by media alone versus in the presence of IL-12 alone with no inhibitors added. Results are expressed as percent inhibition of IFN- $\gamma$  and the IC<sub>50</sub> curve was generated using GraphPad Prism software. IC<sub>50</sub> curves shown were generated from one study and are representative of three

replicate studies. Each data point represents the percent inhibition calculated from the average IFN- $\gamma$  concentration in duplicate wells.

### **Supplemental Table 2: Mouse IL-17RA-Fc Biacore Kinetics**


Kinetics analysis of mIL17RA-Fc was conducted by first amine coupling Protein A (Thermo Scientific, Rockford IL) using standard EDC-NHS chemistry to a Series S CM5 chip (GE Healthcare, Piscataway NJ), followed by capture of mIL17RA-Fc. A titration of 200 to 5.1nM murine IL17A (R&D Systems, Minneapolis MN) was injected across the IL-17RA-Fc surface in HBS-EP+ buffer (GE Healthcare, Piscataway NJ). The surface was regenerated between sample injections using 10mM glycine-HCl pH 2.0 (GE Healthcare, Piscataway NJ). Sensorgrams were double reference subtracted and globally fit to a 1:1 Langmuir model. Binding kinetics of the mouse anti-IL-17 antibody for murine IL-17A/A and IL-17A/F was determined using a Biacore T100 (GE Healthcare, Piscataway NJ). Previous Biacore studies indicated no reactivity of this antibody towards murine IL-17F/F. In brief, biotinylated IL-17A or IL-17A/F, expressed in house using HEK cells, was captured on a streptavidin CAPture chip (GE Healthcare, Piscataway NJ), followed by injection of a 200 to 2.5nM titration of anti-IL-17 antibody in HBS-EP+ buffer. Regeneration of the CAPture reagent surface was done according to the manufacturer's guidelines. Sensorgrams were double reference subtracted and locally fit to a 1:1 Langmuir model.

## Supplemental Figures and Tables:

Supplemental Table 1:

Inhibitor	Replicates	ka (1/Ms)	kd (1/s)	KD (M)
ATI-1221	n=2	8.85E+04	2.17E-05	2.45E-10
ATI-1249	n=6	9.20E+04	2.15E-05	2.34E-10

Supplemental Figure 1:



Supplemental Table 2:

Inhibitor	KD (nM)	KD (nM)	KD (nM)
	IL-17A/A	IL-17A/F	IL-17F/F
mIL-17RA-Fc	0.022	ND	ND
Anti-IL-17 Ab	<0.07	<0.06	No Binding